OMNI DENT[®]
We care. You smile.

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☑️ Instructions for use

OmniCeram évolution 2.0 is a **visible light curing nanohybrid composite** for anterior and posterior restorations. It is based on BIS-GMA-resin and inorganic filler particles of 0.05 – 3 µm. The total filler load is 80 %, the total filler volume is 69 %. The composite material meets the requirements of **DIN EN ISO 4049**.

Indications:

Suitable for all cavities.

Side Effects

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of calcium hydroxide material.

Incompatibility with Other Materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the composite. Neither store the composite material in proximity of eugenol containing products, nor let the composite allow coming into contact with materials containing eugenol.

Preparation of Tooth and Color Matching

Clean the tooth with a fluoride-free polishing paste (e.g. flour of pumice with water) prior to preparation and color matching. Ascertain the tooth shade while teeth are still moist and select the appropriate OmniCeram évolution 2.0 shade. OmniCeram évolution 2.0 is shaded according to VITA®-Shades. It is also available in special shades and bleach shades.

Cavity Preparation

After isolation (best with a rubber dam) prepare the cavity with minimal tooth reduction.

For pulp protection areas close to the pulp should be covered with a thin layer of calcium hydroxide material. Condition and apply a bonding agent according to the manufacturer instructions. The total etch technique with Etching Gel and subsequent application of a light cure bonding agent (e.g. OmniCite 2.0) is recommended. Otherwise use a self etching bonding agent (e.g. OmniCite All'N'One) before application of the OmniCeram évolution 2.0.

Application and Curing of OmniCeram évolution 2.0

Apply OmniCeram évolution 2.0 in the selected shade. OmniCeram évolution 2.0 is easy to model. Place it by using a suitable instrument. Matrix strips may be used.

In case of deep and extended cavities best results are obtained with a build-up of a base liner using a light cure resin modified cement or a light cure floatable composite (e.g. OmniCeram évolution flow 2.0), followed by the application of OmniCeram évolution 2.0.

In case of small cavities OmniCeram évolution 2.0 can be applied directly.

For an optimum result apply OmniCeram évolution 2.0 in **layers of max. 2 mm**. Light cure each increment separately.

By using a polymerization unit (wavelength range 400-500 NM) **with a light intensity of at least 1000mW/cm²**, cure each increment as follows: Lighter shades (e.g. A1, A2, A3, B1, B2, Bleach) **20 seconds** Darker shades (e.g. A3.5, C3, grey) **30 seconds**

Finishing

Depending on the shape of the area contour with finishing diamonds, flexible grinding discs or carbide burs.

Polish to high gloss with polishing discs or polishing brushes.

Storage

Do not store above 25°C (77°F). Protect from direct sunlight. Do not use after expiry date.

Close syringe immediately after use to avoid exposure to light.

PL-Tips are for single use only.

Additional Notes

- Do not use any resin to adjust viscosity of composite restorative material.
- Avoid contact with skin, mucous membrane and eyes.
- Unpolymerized OmniCeram évolution 2.0 may have an irritant effect and can lead to sensitization against methacrylates.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.
- VITA® is a registered trade mark of the VITA-Zahnfabrik, Bad Säckingen, Germany.

Warranty

SanaPro Dental GmbH warrants this product will be free from defects in material and manufacture. SanaPro Dental GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and SanaPro Dental GmbH's sole obligation shall be repair or replacement of the SanaPro Dental GmbH product.

Limitation of Liability

Except where prohibited by law, SanaPro Dental GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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Manufacturer: SanaPro Dental GmbH D-27570 Bremerhaven Germany