Omnisept IMP

Abformdesinfektion

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Impression disinfectant

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GB Instructions for use

Type of application:

Preparation for the disinfection of all dental impression materials.

Notes on application:

Immerse the respective impression into the ready-to-use disinfectant solution, and remove again after the recommended application time. The impression must be fully covered by the solution. Immerse alginates and hydrocolloids for a short time only and then let them dry. Afterwards, rinse the impression with water. Do not mix with cleaning agents or other disinfectants. Depending on the degree of contamination, the working solution may need to be replaced (RKI recommendation: at least once per working day). For professional users only.

Spectrum of efficacy:

Bactericidal (incl. MRSA) as per EN 13727, EN 14561, VAH;
yeastcidal (Candida albicans) as per EN 13624, EN 14562, VAH;
tuberculocidal (M. terrae) in a DGHM (German Society for Hygiene and Microbiology) suspension test, limited virucidal (enveloped viruses, e.g. HBV, HIV and HCV) as per RKI/DVV and EN 17111.

Application times:

Bacteria (incl. MRSA) as per VAH, EN 13727 and EN 14561 (dirty conditions): conc. 1 min.
Candida albicans as per VAH, EN 13624 and EN 14562 (dirty conditions): conc. 1 min.
“Limited virucidal” (enveloped viruses e.g. HBV, HIV, HCV) as per DVV/RKI and EN 17111 (clean conditions): conc. 5 min.
TB (M. terrae in a suspension test, clean conditions): conc. 15 min.

Active ingredients per 100 g:

0.24 g Benzyl-C12-16-alkyldimethylammonium chloride,
0.24 g N,N-didecyl-N,N-dimethyl ammonium chloride,
0.24 g C12-14-alkyl[(ethylphenyl)methyl]dimethyl ammonium chloride, 0.35 g chlorhexidine gluconate

Safety instructions:

H 412 Harmful to aquatic life with long lasting effects. Avoid release to the environment. Dispose of contents/container in accordance with all local, regional, national and international regulations. (EC) 648/2004: Contains < 5% cationic surfactants, fragrances.

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PRISMAN GmbH