Omniwipes wet&easy PLUS

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en:

Type of application: As a medical device: Ready-to-use, alcohol-free product for surface disinfection of non-invasive and invasive medical devices. As a biocidal product: Ready-to-use, alcohol-free wet wipes for disinfecting and cleaning surfaces. The impregnated disinfectant wipes are intended for use in the medical sector, e.g. in (dental) surgeries, hospitals, nursing homes, care and rescue services. Disinfection of medical devices: To disinfect, wipe the surface of the medical device with the soaked cloth. Intended for prophylactic and terminal wipe disinfection of the surfaces of non-critical medical devices such as dental treatment units, cuffs of blood pressure monitors, vein tourniquets, as well as semi-critical medical devices such as invasive ultrasound heads (polycarbonate-free). After terminal disinfection of semi-critical medical devices, rinse them with water (at least drinking water quality) to remove any residues of the disinfectant solution. Not suitable for terminal disinfection of critical medical devices! Disinfection of surfaces: To disinfect and clean small surfaces such as bathtubs, sanitary facilities, etc., wipe the surfaces concerned with the soaked cloth.

Instructions for use: Open the resealable dispenser bag at the top and remove a soaked wipe. Then immediately reseal the bag tightly. For labelling purposes, document the opening and expiray date manually on the enclosed plaque or in the field provided for this purpose before opening the bag for the first time. All surfaces to be disinfected must be completely wetted with the disinfectant solution. Allow the product to take effect according to the specified contact times. Before applying to large areas of sensitive and/or unknown materials such as artificial leather, first test for material compatibility on a small concealed area. Do not use for the treatment of polycarbonate and BMC surfaces! Use is intended for professional users only. The wipes can be used for a maximum of 28 days after the bag has been opened for the first time.

Contact times for disinfection: Bacteria and Candida albicans according to EN 13727 and EN 13624 (dirty conditions): Conc. 1 min. Fungi (A. brasiliensis) according to EN 13624 (clean conditions): Conc. 15 min. Mycobcateria(M. Terrae and M. avium) according to EN 14348 (clean conditions): Conc. 30 min. “Virucidal” (enveloped and non-enveloped viruses according to EN 14476), clean conditions: Conc. 2 min. Contact times for terminal disinfection of semi-critical medical devices: Conc. 30 min. Contact times for the soaked cloths: Bacteria and Candida albicans according to EN 16615 and VAH in dirty conditions: 1 min.

Contents: 80 pieces; 25cm x 25cm, soaked with surface disinfection OmniBiozid PLUS

100g of impregnating liquid contain: 0.48g didecyldimethylammonium chloride; 0.2g N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine.

Only dispose of completely empty packaging in the recycling centre. Dispose of filling material residues as per EAK 07 06 99. Regulation (EC) 648/2004(impregnating liquid): Contains: Fragrances

All serious incidents related to the device shall be reported by the user/patient to the manufacturer and to the competent authority of the State in which the user/patient is established.

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PRISMAN GmbH